
seca 704s/703s

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CONTENTS

1. About this document	25	6. Operation	35
2. Device description	25	6.1 Measure height	35
2.1 Intended use	25	Measuring people less than	
2.2 Description of function	25	130.5 cm tall	36
Scale	25	Measuring people more than	
Measuring rod	26	130.5 cm tall	37
2.3 User qualification	26	6.2 Weighing	38
Assembly	26	7. Hygiene treatment	38
Operation	26	7.1 Cleaning	39
3. Safety precautions	27	7.2 Disinfecting	39
3.1 Safety precautions in these		7.3 Sterilizing	40
instructions for use	27	8. Function check	40
3.2 Basic safety instructions	28	9. Troubleshooting	41
4. Device overview	30	10. Servicing	41
4.1 Deviating/additional controls ..	30	11. Technical data, measuring rod	41
4.2 Identification on the device		11.1 General technical data	41
and on the type plate	31	11.2 Metrological data	42
4.3 Identification on the packaging	32	12. Accessory	42
5. Making the device operational	33	13. Disposal	42
5.1 Deviating/additional scope of		13.1 Disposing of the device	42
supply	33	14. Warranty	42
5.2 Fit column	34	15. Declaration of conformity	43
5.3 Fit measuring rod	34		

1. ABOUT THIS DOCUMENT

This document contains **supplementary** information about the **assembly** and **operation** of the **seca 704s/703s** models. This information relates in particular to the measuring rod included in the scope of supply.

In order to use the **seca 704s/703s** models safely and to be able to make full use of their functions, you should also note the standard instructions for use for the **seca 704/703** models which are enclosed.

Keep both documents safe and ensure that both documents are available at all times.

2. DEVICE DESCRIPTION

2.1 Intended use

The seca medical scale is mainly used in hospitals, medical practices and inpatient care facilities in accordance with national regulations.

In conjunction with the measuring rod included in the scope of supply, the **seca 704s/703s** personal scales serve as measuring stations. The measuring station is for conventional determination of weight and height, establishment of general state of nutrition and assists the physician supervising treatment in making a diagnosis or deciding on a course of treatment.

To make an accurate diagnosis, however, other specific examinations have to be ordered by the physician and their results taken into account, in addition to the determination of weight and height.

The medical measuring rod included in the scope of supply is **not** compatible with other variants of **seca 704/703** models or other seca scales.

2.2 Description of function

Scale The scale functions of the **seca 704s** model are identical to those of the **seca 704** model. You can find information about this in the instructions for use supplied for the **seca 704** model.

The scale functions of the **seca 703s** model are identical to those of the **seca 703** model. You can find information about this in the instructions for use supplied for the **seca 703** model.

Measuring rod The measuring rod for the **seca 704s/703s** models is attached to the scale in such a way that the patient's height can be measured while he or she is standing on the weighing platform.

The measuring rod consists of two telescopic elements:

- lower telescopic element:
measuring range 60 cm -130.5 cm
- upper telescopic element:
measuring range 130.5 cm -200 cm

The scales are on opposite sides of the telescopic elements.

The head slide is latched in the upper telescopic element and can be released for measurements on the lower telescopic element. The measuring flap of the head slide can be folded down.

Use the measuring station only for the purpose quoted in the section entitled "Intended use".

2.3 User qualification

Assembly Devices that are shipped partially assembled may only be mounted by sufficiently qualified persons such as specialist dealers, hospital technicians or seca service technicians.

Operation The device may only be operated by healthcare professionals.

3. SAFETY PRECAUTIONS

3.1 Safety precautions in these instructions for use

**DANGER!**

Used to identify an extremely hazardous situation. If you fail to take note of this information, serious irreversible or fatal injuries will occur.

**WARNING!**

Used to identify an extremely hazardous situation. If you fail to take note of this information, serious irreversible or fatal injuries may result.

**CAUTION!**

Used to identify a hazardous situation. If you fail to take note of this information, minor to moderate injuries may result.

NOTICE!

Used to identify possible incorrect usage of the device. If you fail to take note of this information, you may damage the device, or the measured results may be incorrect.

NOTE

Includes additional information about use of the device.

3.2 Basic safety instructions

- ▶ Follow the instructions for use for the **seca 704/703**.
- ▶ Please observe the information in this user manual.
- ▶ Keep this user manual in a safe place. This user manual is part of the product and must be readily available.

CAUTION!

Patient hazard, damage to device

- ▶ Use the supplied assembly material only.
- ▶ Fold down the measuring slide after each measurement.
- ▶ Have the device serviced regularly as described in the respective section of this document.
- ▶ Technical changes to the device are not permitted. The unit contains no user-serviceable parts. Make sure that maintenance and repair are only carried out by an authorised seca service partner. You can find your local service partner at www.seca.com or send an email to service@seca.com.
- ▶ Make sure you only use genuine seca accessories and spare parts. Otherwise the warranty provided by seca will become null and void.

CAUTION!

Patient hazard

In order to avoid misinterpretations, test results for medical use must be displayed and used in SI units (weight: kilogrammes, length: metres) only. Some devices offer the ability to display test results in other units. This is only an additional function.

- ▶ Use the results exclusively in SI units.
- ▶ The use of measurement results in non-SI units is the sole responsibility of the user.


**CAUTION!****Incorrect measurement due to force shunt**

Operation of the stadiometer affects the weight display of the scale. Weight values which are displayed during a height measurement do not correspond to the patient's actual weight.

- ▶ Ensure that you do not touch the scale whilst weight is being measured.
- ▶ Ensure that the patient does not touch either the stadiometer or the column of the scale whilst weight is being measured.
- ▶ Only read off weight values before or after measuring height.

4. DEVICE OVERVIEW

4.1 Deviating/additional controls

No.	Control	Function
1	Heel stop	<ul style="list-style-type: none"> • Permanently fixed to the weighing platform • Allows patient to be perfectly positioned
2	Column	<ul style="list-style-type: none"> • Prepared for attachment of the measuring rod • Can only be attached in the alignment shown
3	Lower telescopic element	Scale points to the weighing platform Measuring range 1: 60 cm - 130.5 cm
4	Upper telescopic element	Scale points towards display head of scale Measuring range 2: 130.5 cm to 200 cm
5	Measuring flap	<ul style="list-style-type: none"> • Head stop for measuring height • Can be folded down

No.	Control	Function
6	Head slide	<ul style="list-style-type: none"> • Latched in upper telescopic element • Is released for measurements on the lower telescopic element
7	Latch	<ul style="list-style-type: none"> • Latches the head slide for measurements in measuring range 2 • Can be released to push the head slide onto the lower telescopic element (measuring range 1)
8	Read-off edge 1	For measured results in measuring range 1
9	Read-off edge 2	For measured results in measuring range 2

4.2 Identification on the device and on the type plate


Text/symbol	Meaning
	Name and address of manufacturer, date of manufacture
	Model number
	Serial number, consecutive
GAL	Value in m/s^2 (verified models) <ul style="list-style-type: none"> • Gravitational acceleration on earth • Depends on the intended location
	Refer to user manual
	Device complies with EC standards and directives. <ul style="list-style-type: none"> • 0123: Notified body medical products
	Do not dispose of device in household waste

4.3 Identification on the packaging

	Protect from moisture
	Arrows indicate top of product. Transport and store in an upright position.
	Fragile Do not throw or drop.
	Permitted min. and max. temperature for transport and storage
	Permitted min. and max. moisture for transport and storage
	Not sterile
	Do not reuse
	Open packaging here
	Packaging material can be disposed of through recycling programs


5. MAKING THE DEVICE OPERATIONAL

5.1 Deviating/additional scope of supply


No.	Component	Pcs.
a	Set screw 80 mm with thread M5 x 21 mm	2
b	Spacer sleeve	2
c	seca 704s/703s measuring rod	1
d	seca 704s/703s column, prepared for attaching the measuring rod	1
e	seca 704s/703s weighing platform with permanently-fixed heel step	1
	Supplementary instructions for use for seca 704s/703s , no illustration	1
	Instructions for use for seca 704/703 , no illustration	1

5.2 Fit column

1. Place the column on the weighing platform as shown in the adjacent illustration:
 - the display head points away from the platform
 - assembly bores for the measuring rod point towards the weighing platform
2. Perform all the remaining assembly steps for the column as described in the instructions for use for the **seca 704/703**:
 - ▶ Screw column tight
 - ▶ Connect display cable
 - ▶ Connect the power supply cable to the electronics box
3. Set up the scale as described in the instructions for use for the **seca 704/703**.
4. Provide the power supply for the scale as described in the instructions for use for the **seca 704/703**.

5.3 Fit measuring rod

1. Pull out the upper telescopic element as far as the 173 cm mark.
The assembly bores in the lower telescopic element are exposed.
2. Align the measuring rod as shown in the adjacent illustration.
3. Put a set screw in the upper bore of the measuring rod.
4. Align a spacer sleeve as shown in the adjacent illustration.
5. Put the spacer sleeve on the set screw.
6. Screw the measuring rod hand-tight to the column using the upper set screw.
7. Repeat steps 3. to 6. for the lower set screw.
8. Tighten up both set screws firmly.

6. OPERATION


6.1 Measure height

CAUTION!

Bodily injury


- ▶ Perform a function check as described in the section entitled "Function check" on page 40 before using the device each time.

Measuring people less than 130.5 cm tall

1. Release the head slide by pressing on the latch.
2. Fold up the measuring flap until it engages in a horizontal position.
3. Ask the patient to position him or herself under the measuring flap:
 - back to the measuring rod
 - standing up straight
 - heels on the heel stop
 - head straight: Frankfurt Line parallel to the measuring flap
4. Push the head slide on the lower telescopic element down until the measuring flap is in contact with the patient's head.
5. Read the measured result off the read-off edge (read-off edge 1) of the head slide (in this case: 113.5 cm).
6. Ask the patient to step away from the measuring rod.
7. Fold down the measuring flap.
8. Push up the head slide until it engages in the latch.

Measuring people more than 130.5 cm tall

1. Fold up the measuring flap until it engages in a horizontal position.
2. Pull out the upper telescopic element until the patient can get under the measuring flap comfortably.
3. Ask the patient to position him or herself under the measuring flap:
 - back to the measuring rod
 - standing up straight
 - heels on the heel stop
 - head straight: Frankfurt Line parallel to the measuring flap
4. Push the upper telescopic element down until the measuring flap is in contact with the patient's head.
5. Read the measured result off the read-off edge (read-off edge 2) of the lower telescopic element (in this case: 141.5 cm).
6. Ask the patient to step away from the measuring rod.
7. Push the upper telescopic element into the lowest position.
8. Fold down the measuring flap.

6.2 Weighing

CAUTION!

Incorrect measurement due to force shunt

Operation of the stadiometer affects the weight display of the scale. Weight values which are displayed during a height measurement do not correspond to the patient's actual weight.

- ▶ Ensure that you do not touch the scale whilst weight is being measured.
- ▶ Ensure that the patient does not touch either the stadiometer or the column of the scale whilst weight is being measured.
- ▶ Only read off weight values before or after measuring height.

The scale functions of the **seca 704s** model are identical to those of the **seca 704** model. You can find information about this in the instructions for use supplied for the **seca 704** model.

The scale functions of the **seca 703s** model are identical to those of the **seca 703** model. You can find information about this in the instructions for use supplied for the **seca 703** model.

7. HYGIENE TREATMENT

WARNING!

Electric shock

- ▶ Set up the device so that the power supply socket is easy to reach and the device can be disconnected from the power supply quickly.
- ▶ Ensure that your local power supply matches the information on the power supply unit.
- ▶ Do not touch the power supply unit with wet hands.
- ▶ Do not use extension cables or power strips.
- ▶ Make sure that cables are not pinched or damaged by sharp edges.
- ▶ Make sure that cables do not come into contact with hot objects.
- ▶ Do not operate the device at an altitude of more than 3000 m above sea level.

**WARNING!****Electric shock**

The device is not de-energized when the on/off key is pressed and the display goes out. Use of fluids on the device may cause electric shock.

- ▶ Ensure that the device is switched off before performing any hygiene treatment.
- ▶ Disconnect the power supply connector before performing any hygiene treatment.
- ▶ Before each hygiene treatment, take the rechargeable battery out of the device (if present and removable).
- ▶ Ensure that no fluids penetrate the device.

**CAUTION!****Damage to device**

Inappropriate detergents and disinfectants may damage the sensitive surfaces of the device.

- ▶ Use only disinfectants suitable for sensitive surfaces (e.g. active ingredients including 70 % isopropanol).
- ▶ Do not use aggressive or abrasive cleaning agents.
- ▶ Do not use organic solvents (e.g. white spirit or petroleum spirit).

7.1 Cleaning

- ▶ Use a soft cloth dampened with mild soapsuds to clean the surfaces of the device.

7.2 Disinfecting

1. Check that your disinfectant is suitable for sensitive surfaces and acrylic sheet.
2. Follow the instructions for use for the disinfectant.
3. Disinfect the device:
 - ▶ moisten a soft cloth with disinfectant and wipe down the device with it.
 - ▶ comply with the intervals: see table.

Interval	Component
Before Each measurement	<ul style="list-style-type: none"> • Weighing platform • Head slide with measuring flap
After Each measurement	<ul style="list-style-type: none"> • Weighing platform • Head slide with measuring flap
As required	<ul style="list-style-type: none"> • Column and display head of scale • Telescopic elements of the stadiometer • Measuring flap extension (if present)

7.3 Sterilizing

This device may not be sterilised.

8. FUNCTION CHECK

- ▶ Perform a function check prior to each use.

A complete function check includes:

- visual inspection for mechanical damage
- checking the alignment of the device
- visual and function check of the display elements
- function check of all the controls shown in the section entitled “Overview”
- function check of optional accessories

If you notice any faults or deviations during the function check, first try to resolve the error with the aid of the section entitled “What do I do if ...” in this document.

CAUTION! **Personal injury**

If you notice any faults or deviations during the function check which cannot be resolved with the aid of the section entitled “What do I do if ...” in this document, you may not use the device.

- ▶ Have the device repaired by seca service or by an authorized service partner.
- ▶ Follow the section entitled “Servicing” in this document.

9. TROUBLESHOOTING

Problem	Cause/solution
... the measurement is not a plausible height?	Measured result read off incorrectly <ul style="list-style-type: none"> - Use the correct read-off edge for the current measuring range - Read off decimal places for the current measuring range in the correct direction Patient in wrong position/posture <ul style="list-style-type: none"> - Correct patient's position/posture

10.SERVICING

The seca measuring rod does not require any servicing. However, if anything ever fails to work, contact an authorized service partner.

You can find a service partner in your area at www.seca.com or by sending an e-mail to service@seca.com

11.TECHNICAL DATA, MEASURING ROD

11.1 General technical data

General technical data	
Dimensions <ul style="list-style-type: none"> • Depth • Width • Height (retracted) • Height (extended) 	37 mm 57 mm 1324 mm 2320 mm
Net weight	approx. 1 kg
Ambient conditions, operation <ul style="list-style-type: none"> • Temperature • Air pressure • Humidity 	+10 °C to +40 °C (50 °F to 104 °F) 700 hPa - 1060 hPa 30 % - 80 %, no condensation
Ambient conditions, storage <ul style="list-style-type: none"> • Temperature • Air pressure • Humidity 	-10 °C to +65 °C (14 °F to 149 °F) 700 hPa - 1060 hPa 0 % - 95 %, no condensation
Ambient conditions, transport <ul style="list-style-type: none"> • Temperature • Air pressure • Humidity 	-10 °C to +65 °C (14 °F to 149 °F) 700 hPa - 1060 hPa 0 % - 95 %, no condensation

General technical data	
Medical device in accordance with Directive 93/42/EEC	Class I with measuring function

11.2 Metrological data

Metrological data	
Measuring range 1	60 cm - 130.5 cm/ 23.6 inch - 51.4 inch
Measuring range 2	130.5 cm - 200 cm/ 51.4 inch - 78.7 inch
Graduations	1 mm/1/8 inch
Accuracy	± 5 mm/± 0.2 inch

12.ACCESSORY

seca 430 measuring flap	430-00-00-009
-------------------------	---------------

13.DISPOSAL

13.1 Disposing of the device

Do not dispose of the device with household waste. The device must be disposed of properly as electronic waste. Comply with the national provisions applicable in your country. For further information contact our service department at:

service@seca.com

14.WARRANTY

We offer a two-year warranty from the date of delivery for defects attributable to faulty material or poor workmanship. This excludes all moveable parts such as (rechargeable) batteries, cables, power supply units, etc. Defects which are covered by the warranty shall be rectified free of charge for customers on production of the sales receipt. No further claims can be accepted. The costs of shipment in both directions shall be borne by the customer where the device is not located at the customer's premises. In the event of any damage during shipment warranty claims can only be asserted where the complete original packaging was used for

shipment and the scales were secured inside in the same manner as in the original packaging. You should therefore keep all packaging.

The warranty shall become null and void where the device is opened by persons not expressly authorised to do so by seca.

In the event of a warranty issue, please contact your local seca office or the dealer from whom you ordered the product.

15. DECLARATION OF CONFORMITY

seca gmbh & co. kg hereby declares that the product meets the terms of the applicable European directives. The unabridged declaration of conformity can be found at: www.seca.com.

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