
seca 899

Instructions for use

17-10-07-576-100e_2025-02S

Deutsch – Gebrauchsanweisung	3
English – Instructions for use	33
Français – Mode d’emploi.....	62
Italiano – Istruzioni per l’uso	94
Español – Instrucciones de uso	124
Português – Manual de instruções de utilização	156
Polski – Instrukcja obsługi.....	187
Dansk – Brugsanvisning.....	218
Svenska – Bruksanvisning.....	246
Norsk – Bruksanvisning.....	274
Suomi – Käyttöohje	302
Nederlands – Gebruiksaanwijzing	330
Ελληνικά – Οδηγίες χρήσης.....	360

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CONTENTS

1	Description of device	34		
1.1	Intended use.....	34		
1.2	User qualification.....	34		
1.3	Patient target group.....	34		
1.4	Contraindications	34		
1.5	Clinical benefit.....	34		
1.6	Functional description	34		
2	Safety information	35		
2.1	Safety information in these instructions for use .	35		
2.2	Basic safety information ...	35		
	Handling device.....	35		
	Preventing electric shock.....	36		
	Preventing injuries and infections	37		
	Preventing damage to device	38		
	Handling measuring results	39		
	Handling packaging material.....	39		
	Handling batteries and rechargeable batteries..	40		
3	Overview	41		
3.1	Controls.....	41		
3.2	Symbols in the display	42		
3.3	Labels	43		
4	Starting up device.....	46		
4.1	Scope of delivery.....	46		
4.2	Table stand.....	47		
4.3	Installing the wall bracket .	47		
4.4	Establishing power supply	48		
	Inserting batteries.....	48		
	Connecting the power supply unit (optional)	48		
4.5	Connecting the connecting cable to the weighing platform	49		
4.6	Setting up and aligning the device	49		
5	Operation.....	50		
5.1	Switching the device on and off.....	50		
	Switching on the device	50		
	Switching off the device	50		
5.2	Weighing.....	51		
	Weighing a patient.....	51		
	Determining Body Mass Index (BMI).....	51		
	Permanently displaying the weight (Hold).....	52		
	Taring additional weight (Tare).....	52		
6	Hygiene treatment	54		
6.1	Cleaning.....	54		
6.2	Disinfecting.....	55		
6.3	Sterilizing	55		
7	Function check.....	55		
8	Troubleshooting.....	56		
9	Servicing	58		
10	Technical data.....	58		
10.1	General technical data.....	58		
10.2	Weight measurement.....	59		
11	Optional accessories.....	60		
12	Disposal	60		
12.1	Disposing of the device	60		
12.2	Disposing of batteries and rechargeable batteries.....	60		
13	Warranty	61		
14	Declaration of conformity	61		

1 DESCRIPTION OF DEVICE

1.1 Intended use

The electronic flat scale supports physicians in decision-making regarding diagnosis or therapy, based on the patient's weight or weight changes.

1.2 User qualification

The device may only be used by persons with sufficient expertise.

1.3 Patient target group

The scale is intended for persons of any age who can stand independently on the scale and who do not exceed the maximum weight capacity of the scale.

1.4 Contraindications

No contraindications are known.

1.5 Clinical benefit

The electronic flat scale supports physicians in decision-making regarding diagnosis or therapy, based on the measured weight.

By that, the electronic flat scale only provides an indirect clinical benefit.

The physician must perform additional examinations for an accurate diagnosis.

1.6 Functional description

Weight on this electronic personal scale is detected by four load cells.

Body Mass Index (BMI) can be calculated from weight and height. Height is entered manually.

2 SAFETY INFORMATION

2.1 Safety information in these instructions for use

**DANGER!**

Used to identify an extremely hazardous situation. If you fail to take note of this information, serious irreversible or fatal injuries will occur.

**WARNING!**

Used to identify an extremely hazardous situation. If you fail to take note of this information, serious irreversible or fatal injuries may result.

**CAUTION!**

Used to identify a hazardous situation. If you fail to take note of this information, minor to moderate injuries may result.

NOTICE!

Used to identify possible incorrect usage of the device. If you fail to take note of this information, you may damage the device, or the measuring results may be incorrect.

NOTE

Includes additional information about use of the device.

2.2 Basic safety information

Handling device

- ▶ Please take note of the information in these instructions for use.
- ▶ Keep the instructions for use in a safe place. The instructions for use are a component of the device and must be available at all times.
- ▶ In the interest of patient safety, you and your patients are obliged to report serious events that occur in connection with this product to the manufacturer and to the authority responsible in your country.

DANGER!
Risk of explosion

- ▶ Do not use the device in an environment in which one of the following gases has accumulated:
 - oxygen
 - flammable anesthetics
 - other flammable substances/air mixtures

CAUTION!
Patient hazard, damage to device

- ▶ Have servicing carried out regularly as described in the relevant section of this document.
- ▶ Technical modifications may not be made to the device. The device does not contain any user-serviceable parts. Only have servicing and repairs performed by an authorized seca service partner. You can find your local service partner at www.seca.com.
- ▶ Use only seca original accessories and spare parts, otherwise seca will not grant any warranty.

CAUTION!
Patient hazard, malfunction

- ▶ Keep medical electrical devices such as high-frequency surgical equipment at a minimum distance of approx. 1 meter to prevent faulty measurements.
- ▶ Keep high-frequency equipment such as cell-phones at a minimum distance of approx. 1 meter to prevent faulty measurements.
- ▶ The actual transmission output of HF devices may require minimum distances of more than 1 meter. For details, go to www.seca.com.

Preventing electric shock

WARNING!
Electric shock

- ▶ Set up devices which can be operated with a power supply unit so that the power supply socket is within easy reach and the power supply can be disconnected quickly.
- ▶ Ensure that your local electricity supply matches the details on the power supply unit.

- ▶ Never touch the power supply unit with wet hands.
- ▶ Do not use extension cables and multiple outlets.
- ▶ Ensure that cables are not crushed or damaged by sharp edges.
- ▶ Ensure that cables do not come into contact with hot objects.
- ▶ Do not operate the device at an altitude of more than 3000 m above sea level.

Preventing injuries and infections

WARNING! **Injury from falling**

- ▶ Ensure that the device is steady and level.
- ▶ Route connector cables (if present) so that neither users nor the patient can trip over them.
- ▶ Assist people with limited motor skills when they are getting up, e.g. from a wheelchair.
- ▶ Ensure that the patient does not step directly onto or off the edges of the weighing platform.
- ▶ Ensure that the patient steps onto and off the weighing platform slowly and safely.

WARNING! **Danger of slipping**

- ▶ Ensure that the weighing platform is dry before the patient steps onto it.
- ▶ Ensure that the patient has dry feet before stepping onto the weighing platform.
- ▶ Ensure that the patient steps onto and off the weighing platform slowly and safely.

WARNING! **Risk of infection**

- ▶ Before and after every measurement, wash your hands to reduce the risk of cross-contamination and nosocomial infections.
- ▶ Should the patient have contagious diseases, subject the device to a hygiene treatment as described in the relevant section of these instructions for use immediately afterwards.

Preventing damage to device

- ▶ Ensure that the patient does not have any open wounds or infectious skin alterations which may come into contact with the device.
- ▶ Subject the device to a hygiene treatment at regular intervals.

NOTICE!

Damage to device

- ▶ Ensure that fluids never get inside the device. These can destroy the electronics.
- ▶ For devices with power supply operation: Switch off the device before you disconnect the power supply unit from the power supply socket.
- ▶ For devices with power supply operation: If the device is not to be used for an extended period, disconnect the power supply unit from the power supply socket. Only then is the device de-energized.
- ▶ For devices with battery or rechargeable battery operation: If you are not using the device for an extended period of time, remove batteries or rechargeable batteries. Only then is the device de-energized.
- ▶ Do not drop the device.
- ▶ Do not subject the device to shocks or vibrations.
- ▶ Perform a function check before each use as described in the corresponding section in this document. Do not operate the device if it is not working properly or is damaged.
- ▶ Do not place the device in direct sunlight and ensure that it is not placed in the direct proximity of a heat source. The excessive temperatures could damage the electronics.
- ▶ Avoid rapid temperature fluctuations. If the device is transported so that a temperature difference of over 20 °C occurs, the device must be left to stand for at least 2 hours before it is switched on, otherwise condensation may form; this may damage the electronics.

- ▶ Use the device only in the intended ambient conditions outlined in the section entitled “Technical data”.
- ▶ Store the device only in the intended storage conditions outlined in the section entitled “Technical data”.
- ▶ Use only cleaning agents and disinfectants which match the details in the section entitled “Hygiene treatment”.
- ▶ For scales: Ensure that maximum capacity is not exceeded.

Handling measuring results

NOTICE!

Inconsistent measuring results

- ▶ Before you save and continue using measured values determined using this device (e.g. in seca software or in an EMR system), ensure that the measured values are plausible.
- ▶ If measured values have been sent to seca software or to an EMR system, ensure before continuing to use them that the measured values are plausible and assigned to the correct patient.

Handling packaging material

WARNING!

Risk of suffocation

Packaging materials made of plastic film (bags) present a risk of suffocation.

- ▶ Store packaging material out of the reach of children.
- ▶ If the original packaging material is no longer available, only use plastic bags with safety holes to reduce the risk of suffocation. Use recyclable materials if possible.

NOTE

Store the original packaging material for future use (e.g. returning for servicing).

Handling batteries and rechargeable batteries

WARNING!

Injury

Batteries and rechargeable batteries contain harmful substances which may explode and so escape if not handled properly.

- ▶ Do not try to recharge batteries.
- ▶ Do not expose batteries/rechargeable batteries to heat.
- ▶ Do not burn batteries/rechargeable batteries.
- ▶ If harmful substances are leaking out, avoid contact with the skin, eyes, and mucous membranes. Rinse affected areas of the body with plenty of clean water and seek medical help at once.


NOTICE!

Damage to device and malfunction due to incorrect handling

- ▶ Use only the type of battery/rechargeable battery quoted in this document.
- ▶ Always replace all batteries/rechargeable batteries at the same time.
- ▶ Do not short-circuit batteries/rechargeable batteries.
- ▶ If you do not use the device for an extended period of time, remove the batteries/rechargeable batteries. This prevents acid from leaking into the device.
- ▶ If harmful substances have penetrated the device, do not continue using it. Have the device checked by an authorized seca service partner and repaired if necessary.


3 OVERVIEW

3.1 Controls

Item	Control	Function
1		Switches scale on and off
2		<ul style="list-style-type: none"> • Sets the display to zero • Activates/deactivates Tare function • In the menu: Increases value
3		<ul style="list-style-type: none"> • Activates/deactivates Hold function • In the menu: Reduces value
4		<ul style="list-style-type: none"> • Activates BMI function • In the menu: Confirms value
5	Foot screw	4 pcs, for precise alignment
6	Power supply connection	For connecting a power supply unit (available as an option)
7	Battery compartment	Takes 6 AA-type batteries (AA/1.5 V)
8	Spirit level	Indicates whether the device is horizontal




3.2 Symbols in the display

Item	Symbol	Meaning
1	kg	Weight value in kilograms
2	NET	Tare function active
3	HOLD	Hold function active
4		Non-verifiable function active
5	BMI	BMI function active
6		Power supply unit connected
7		Batteries are low

3.3 Labels


Markings on the device and on the type plate	
Text/symbol	Meaning
	Name and address of manufacturer, date of manufacture
UDI	Unique Device Identification (product identification number as per Directive (EU) 2017/745)
	Article number
	Serial number
GAL	Value in m/s^2 (model-dependent) <ul style="list-style-type: none"> • Gravitational acceleration on earth • Depends on the intended location
ProdID	seca product identification number
Approval Type	Type designation of design approval
Mat.No.	Variant number
	Follow instructions for use
	Medical electrical device, Type B
IP20	Type of protection to IEC 60529: <ul style="list-style-type: none"> • Protection against ingress of solid foreign bodies with a diameter of over 12.5 mm • Protection against access with fingers • No protection from ingress of water
e	Value in units of mass Used to classify and verify a scale
	Class III scale in accordance with Directive 2014/31/EU

Markings on the device and on the type plate	
Text/symbol	Meaning
	<p>Device complies with EU directives</p> <ul style="list-style-type: none"> • M: Conformity label in compliance with directive 2014/31/EU for non-automatic scales (verified models) • 24: (Example: 2024) Year in which declaration of conformity was implemented and the CE label was applied (verified models) • 0102: Notified Body for Metrology (verified models) • 0123: Notified Body for Medical Devices
	<p>Medical device in accordance with Regulation (EU) 2017/745</p>
	<p>Device complies with United Kingdom directives</p> <ul style="list-style-type: none"> • M: Conformity label in compliance with UK directive SI 2016 no. 1152 for non-automatic scales (NAWIR) (verified models) • 24: (Example: 2024) Year in which declaration of conformity was implemented and the UKCA label was applied (verified models) • xxxx: Notified Body for Medical Devices of the United Kingdom • yyyy: Notified Body for Metrology of the United Kingdom (verified models)
 	<p>Importer/representative in the United Kingdom: seca Ltd 40 Barn Street B5 5QB Birmingham United Kingdom</p>
 	<p>Importer/representative in Switzerland: seca ag (schweiz) Medizinische Waagen und Messsysteme Schönmatt Str. 2 CH-4153 REINACH</p>
	<p>Type plate on the power supply connection socket</p> <ul style="list-style-type: none"> • Power supply voltage required in V • Maximum current consumption in mA • : Note polarity of device connector • : Operate device with direct current • : Use compatible seca power supply units only
	<p>Do not dispose of device in household waste</p>

Markings on the packaging	
Symbol	Meaning
	Protect from moisture
	Arrows indicate top of product Transport and store in an upright position
	Fragile, do not throw or drop
	Permitted min. and max. temperature for transport and storage
	Permitted min. and max. humidity for transport and storage
	Permitted min. and max. air pressure for transport and storage
	Open packaging here
	Packaging material can be disposed of through recycling programs


4 STARTING UP DEVICE

4.1 Scope of delivery

Item	Component	Pcs.
a	Weighing platform	1
b	Display unit with connecting cable	1
c	Wall plug	2
d	Wall bracket	1
e	Screw	2
f	Battery	6
-	Instructions for use	1


4.2 Table stand


The display unit can be put on a stand rather than mounted on a wall.

- ▶ To make the display as easy as possible to read, swivel the display unit to the desired position.


4.3 Installing the wall bracket

1. Screw the wall bracket to the wall as shown. Note that the cable between the scale and the display unit is 2 m long.


2. Attach the display unit to the wall bracket.

3. To transport the scale, remove the display unit from the wall bracket and push it into the holder under the weighing platform.

4.4 Establishing power supply

Inserting batteries

You need 6 batteries (type AA/1.5 V, included in scope of delivery).

1. Press on the latch of the battery compartment.
2. Remove the cover of the battery compartment.

NOTICE!

Damage to device

Device may overheat and become damaged.

- ▶ Ensure the correct polarity of the batteries (markings in the battery holder).

3. Insert the batteries.

NOTE

If **bAtt** appears in the display, you have inserted one of the batteries the wrong way round or the batteries have discharged.

4. Close the battery compartment.

Connecting the power supply unit (optional)

The scale can be operated with a power supply unit (optional accessory).


WARNING!

Personal injury or damage to the device as a result of incorrect power supply units


Conventional power supply units may deliver a higher voltage than is indicated on them. The measuring device may overheat, catch fire, melt or short-circuit.

- ▶ Use only original seca power supply units with a controlled 12 V output voltage.

1. Insert the connector required for your power supply into the power supply unit.
2. Insert the device connector of the power supply unit into the connecting socket of the scale.
3. Plug the power supply unit into a power supply socket.

4.5 Connecting the connecting cable to the weighing platform

1. Connect the cable to the connection underneath the carrying handle.
2. Set the scale upright.

NOTICE!

Faulty measurement

If the display cable is touching the weighing platform, the patient's weight cannot be measured correctly.

- ▶ Route the cable so that it cannot touch the weighing platform.

4.6 Setting up and aligning the device


CAUTION!

Faulty measurement as a result of devices being set up incorrectly

Soft floors (wooden boards, for example) give under the patient's weight and falsify the measuring result. Incorrect ambient conditions or lack of alignment lead to measurement errors.

- ▶ Select a setup location with a level, stable floor in order to achieve precise measuring results.
- ▶ Use the device only in the intended ambient conditions outlined in the section entitled "Technical data".
- ▶ Align the device before using it and every time it changes location.

1. Place the device on a firm, level surface.
2. Align the device by turning the foot screws so that the air bubble of the spirit level is in the exact center of the circle.
3. Ensure that only the feet of the device are in contact with the ground.

5 OPERATION

CAUTION!

Personal injury/faulty measurement

- ▶ Before using the device each time, perform a function check as described in the corresponding section of these instructions for use.

5.1 Switching the device on and off

Switching on the device

- ✓ The scale has no load.
- ▶ Press the **start** key.
 - ⇒ **SECA** appears in the display.
 - ⇒ All the elements of the display are shown briefly.
 - ⇒ The scale is operational when **0.0** appears in the display.
 - ⇒ If the scale is operated with a power supply unit, the  symbol appears in the display.

Switching off the device

- ▶ Press the **start** key.
- ⇒ The display goes out.

NOTE

For battery pack/battery mode only: The device switches off automatically after a short time if there is no load on the scale.

5.2 Weighing

Weighing a patient

- ✓ The scale has no load.
- 1. If the scale switched off automatically after the last measurement, press the **start** key.

NOTICE!

Damage to device

The display is sensitive and may be damaged under load

- ▶ Make sure that the patient does not step on the display.

- 2. Ask the patient to step onto the scale.
- 3. Read off the measuring result.

Determining Body Mass Index (BMI)

Body Mass Index correlates height and weight.

- ✓ The scale has no load.

- 1. Briefly press the **bmi** key.
 - ⇒ The last height entered appears on the display.
- 2. Adopt the saved value (**bmi** key) or change the value using the arrow keys and then confirm with the **bmi** key.
 - ⇒ The **BMI** function is activated.
- 3. Ask the patient to step onto the scale and keep still.
 - ⇒ The patient's BMI is automatically calculated and displayed.
- 4. Read off the BMI.
- 5. To deactivate the **BMI** function, press the **bmi** key.

Permanently displaying the weight (Hold)

If the **Hold** function is activated, the weight continues to be displayed once the load is removed from the scale. This enables you to attend to the patient before recording the weight value.

✓ The scale has no load.

1. Ask the patient to step onto the scale.

2. Press the **hold** key.

3. Wait until the weight value has stopped flashing.

⇒ The  symbol and the **HOLD** message are displayed.

4. To deactivate the function, press the **hold** key.

⇒ The  symbol and the **HOLD** message are no longer displayed.

Taring additional weight (Tare)

Use the **Tare** function to prevent an additional object (e.g. walking stick) affecting the patient's weight value.

WARNING!

Faulty measurement with the Tare function

If the **Tare** function is used incorrectly, the weight of the patient may not be determined correctly.

► Perform the steps as described to prevent faulty measurements.

CAUTION!

Faulty measurement as a result of force shunt

If an additional weight is touching the surface on which the device is standing, the weight will not be measured correctly.

► Make sure that additional weights are placed only on the weighing platform.

NOTE

The maximum weight which can be displayed is reduced by the weight of the objects which have been tared.

✓ The scale has no load.

1. Place an object on the weighing platform.
⇒ The object's weight is displayed.
2. Press the **tare** key.
⇒ The display is reset to zero.
⇒ **0.0** flashes in the display until the weight is detected to be stable.
⇒ The **NET** message is displayed.

NOTE

If you remove the additional weight again, - - - - - is displayed.

3. Ask the patient to step onto the weighing platform with the additional weight (or carrying the additional weight).
⇒ The patient's weight is displayed. The additional weight is deducted automatically.
4. Read off the measuring result.
5. To deactivate the **Tare** function, press the **tare** key.
⇒ The **NET** message is no longer displayed.

NOTE

If the device is switched off, the **Tare** function is deactivated automatically.

6 HYGIENE TREATMENT

WARNING! **Electric shock**

The device is not de-energized when the on/off button is pressed and the display goes out. Use of fluids on the device may cause an electric shock.

- ▶ Before each hygiene treatment, ensure that the device is switched off.
- ▶ Disconnect the power supply connector before each hygiene treatment.
- ▶ Take the batteries out of the device before every hygiene treatment.
- ▶ Ensure that no fluids penetrate the device.

WARNING! **Risk of infection**

- ▶ Subject the device to a hygiene treatment at regular intervals as described in this section.

NOTICE!

Damage to device

Unsuitable cleaning agents and disinfectants may damage the sensitive surfaces of the device and lead to transparent components becoming clouded.

- ▶ Do not use aggressive or abrasive cleaning agents.
- ▶ Do not use organic solvents (e.g. white spirit or petroleum spirit).

6.1 Cleaning

- ▶ If required, moisten a soft cloth with a mild soap solution and wipe the device over with it.

6.2 Disinfecting

1. Disinfect the device at regular intervals using a disinfectant suitable for sensitive surfaces and acrylic glass. seca recommends a 70 % ethanol solution as a surface disinfectant.
2. Follow the instructions for use of the disinfectant.
3. Disinfect the device as follows:
 - ▶ Moisten a soft cloth with disinfectant and wipe down the device with it.
 - ▶ Comply with the intervals, see table:

Interval	Component
Before and after every measurement with direct skin contact	Weighing platform
As required	Display unit

6.3 Sterilizing

The device must not be sterilized.

7 FUNCTION CHECK

- ▶ Perform a function check before each use.

A complete function check includes:

- A visual inspection for mechanical damage
- A test of the alignment of the device
- A visual and function check of the display elements
- A function check of all the controls shown in the section entitled “Overview”
- A function check of the optional accessories

If you find faults or deviations during the function check, first try to remedy the fault with the aid of the “Troubleshooting” section in this document.

CAUTION!
Personal injury

If you find faults or deviations during the function check which you are unable to remedy with the aid of the “Troubleshooting” section in this document, you must not use the device.

- ▶ Have the device repaired by seca Service or an authorized service partner.
- ▶ Follow the section entitled “Servicing” in this document.

8 TROUBLESHOOTING

If faults occur when operating the device, first attempt to remedy them yourself using the following table. If the fault persists, contact seca Service or an authorized service partner.

Fault	Cause	Remedy
The scale does not display a weight although a load was placed on it	Scale is switched off or without batteries or not correctly connected	<ul style="list-style-type: none"> • Check that the scale is switched on • Check that batteries are inserted • Check that the cable for the display unit is correctly connected • Check that the power supply unit is correctly connected
0.0 does not appear before weighing	Scale had load on it before being switched on	<ul style="list-style-type: none"> • Remove the load from the scale • Switch off the scale, then switch it back on again
Single segment is illuminated continuously or not at all	Display faulty	Inform seca Service
Display: 	Battery voltage is dropping	Replace the batteries
Display: bAtt	Batteries discharged or not inserted	<ul style="list-style-type: none"> • Replace the batteries • Insert the batteries the right way round

Fault	Cause	Remedy
Display: StOP	Maximum capacity exceeded	Remove the load from the scale
Display flashing	Consistent weight value not yet detected	<ul style="list-style-type: none"> • One of the following functions is active: Hold, Autohold, Tare. Wait until the device has detected a consistent weight value. The display stops flashing. • Take the load off the device, wait until 0.0 is displayed and repeat the measurement procedure.
Display: E + figure	Too high a load on the scale or too high a load in one corner	<ul style="list-style-type: none"> • Remove the load from the scale or distribute the weight more evenly • Switch off the scale, then switch it back on again • Remove batteries and put back in • If the problem persists, inform seca Service

9 SERVICING

CAUTION!

Faulty measurements as a result of verification being performed poorly or not at all

- ▶ Have verification performed only by authorized persons.
- ▶ Always have verification performed if one or more quality seals are damaged.

seca recommends having your device serviced prior to verification.

CAUTION!

Faulty measurements as a result of poor servicing

- ▶ Have servicing and repairs carried out exclusively by seca Service or an authorized service partner.
- ▶ You can find your local service partner at www.seca.com.

10 TECHNICAL DATA

10.1 General technical data

Dimensions of weighing platform <ul style="list-style-type: none">• Depth• Width• Height	356 mm 321 mm 59 mm
Dimensions of display unit <ul style="list-style-type: none">• Depth• Width• Height	47 mm 205 mm 73 mm
Net weight	Approx. 4.3 kg/ 9.5 lbs
Ambient conditions, operation <ul style="list-style-type: none">• Temperature• Air pressure• Humidity	+10 °C to +40 °C (50 °F to 104 °F) 700 – 1060 hPa 30 % – 80 %, no condensation

17-10-07-576-100e_2025-02S

Ambient conditions, storage/transport <ul style="list-style-type: none"> • Temperature • Air pressure • Humidity 	-10 °C to +65 °C (14 °F to 149 °F) 700 – 1060 hPa 0 % – 95 %, no condensation
Height of digits	20 mm
Power supply <ul style="list-style-type: none"> • Batteries <ul style="list-style-type: none"> – Supply voltage – Battery type • Power supply unit (optional) <ul style="list-style-type: none"> – Supply voltage – Insulated device, protection class II (EN 60601-1) 	9 V 6 x type AA, 1.5 V, not rechargeable 12 V 
Maximum current consumption	Typically 30 mA
Medical device in accordance with Regulation (EU) 2017/745	Class I with measuring function
EN 60601-1: Medical electrical equipment, type B	
Application part in line with IEC 60601-1	Tread area
Type of protection	IP20
Duty cycle	Continuous duty

10.2 Weight measurement

Verification in line with Directive 2014/31/EU	Class III
Maximum capacity	200 kg
Accuracy <ul style="list-style-type: none"> • 0 to 50 kg • 50 kg to 150 kg • 150 kg to 200 kg 	± 50 g ± 100 g ± 200 g
Graduation <ul style="list-style-type: none"> • 0 to 150 kg • 150 kg to 200 kg 	100 g 200 g

11 OPTIONAL ACCESSORIES

Accessory	Article number
seca 421 carry case	421 00 00 009
Power supply: Switch-mode power supply unit: 100-240 V~/50-60 Hz/12 V=/min. 0.5 A	68 32 10 265 009 or 68 32 10 270 009

12 DISPOSAL

12.1 Disposing of the device

Do not dispose of the device in your household waste. The device must be properly disposed of as electronic scrap. Follow your respective national regulations. For more information, please contact seca Service or an authorized service partner.

12.2 Disposing of batteries and rechargeable batteries

Do not dispose of used batteries and rechargeable batteries in household waste, regardless of whether they contain harmful substances or not. As a consumer, you are legally obliged to dispose of batteries and rechargeable batteries via local authority collection points or trade collection points. Only dispose of batteries and rechargeable batteries once they are completely discharged.

13 WARRANTY

There is a two-year warranty period from delivery for defects attributable to poor materials or workmanship. All movable parts, e.g. batteries, cables, power supply units, rechargeable batteries etc. are exempt. Defects which come under the warranty will be repaired for the customer free of charge against proof of purchase. Additional claims cannot be considered. Costs of transport to and from seca are the responsibility of the customer if the device is located somewhere other than the customer's headquarters. In the event of transport damage, claims under warranty can only be made if the complete original packaging was used for transport and the device was secured and fastened in it according to its originally packaged condition. You should therefore keep all packaging parts.

The warranty will be voided if the device is opened by persons not expressly authorized by seca to do so.

In the event of a warranty issue, please contact your local seca office or the dealer from whom you ordered the product.

14 DECLARATION OF CONFORMITY

seca gmbh & co. kg hereby declares that the product complies with the terms of the applicable European directives and regulations. The unabridged declaration of conformity can be found at www.seca.com.

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